

BD FACSAriaFusion™ basic SOP

KEY SPECIFICATIONS

Low to high sample pressure (10 min)

Amplitude: intensity of the drive frequency: drops formed per second

Drop 1: pixels from the top of the first dispersed drop of pixels at the breakoff

Data displayed on 262,144 scale

Fluid cart:

5L tanks

BD FACSAriaFusion™ basic SOP

SYSTEM STARTUP

1.	Check fluids, empty waste, fill sheath and attach sheath
2.	Check the waste container back pressure problems
3.	Turn on BSC and compressor.
4.	Turn on computer and cylinder order.
5.	Start FACSDiva software, select Use Current CST settings
6.	Cytometer > Fluidics startup if a Fluidics Shutdown was otherwise perform Flow Cell with DI.
7.	Clean the plates, stream camera, breakoff camera, the sensor behind the nozzle.

BD FACSAriaFusion™ basic SOP

SYSTEM SHUTDOWN

- | | |
|----|--|
| 1. | Stop stream, remove nozzle, clean and sonicate in DI, stop flow, and install closed loop nozzle. |
|----|--|

cells).

BD FACSAriaFusion™ basic SOP

EXPERIMENT SETUP

- Large or new cells, changes in pressure, particles require different area scaling, area scaling should be performed on cells (large cells have larger area). Draw Area vs. Height plot and draw FSC and each FL. Change on Cytometer window > Laser tab. Use settings to preserve ASF and LD.
- Area >= Height, important for sorting!
- **Sheath pressure > velocity > pulse width > area scaling > voltages > compensation.**
-

BD FACSAriaFusion™ basic SOP

SORTING

- Enrich with a high threshold > re-sort with a lower threshold.
- Donot use exponential display or snap-to gates. All events below 0 in 10^{exp} in a sort gate get sorted into the gate, even those outside the gate.
- Rare events > use storage gate > sort to far right stream > less loss
- Yield mask = drop sorted, reduce to 1 drop sort. for fanning
- Purity mask = incidence
- Phase mask = cell positioning, use to correct for fanning, yields
-

BD FACSAriaFusion™ basic SOP

Fusion sort collection camera and stream targeting

1. Open the StCamSWare application for camera view. Use camera or upper position with the light deflector to tube to collection.
2. Open the IconSortDeviceControl application to operate the Sort Alignment Software.
3. Change the sort device as required.

ACDU

1. Install the splash guard.
2. Clear the sort chamber, and eject the ACDU stage.
3. Choose the sort device or create a custom device. Adjust the height and apply.
4. Sort with the default (50-87s)-8.4 (p17)6843y(a1)E4tc3Q(d)W8

BD FACSAriaFusion™ basic SOP

FAQ & Troubleshooting

- Most problems will be fluidics (air waste/orange, sheath/blue).
- Bubbles, turn stream on and off, perform Fluidics Startup.
- Sort > backflush.
-

BD FACSAriaFusion™ basic SOP

SYSTEM MAINTENANCE

(performed by Flow Cyt Lab technician or designated alternate)

At least TWICE WEEKLY maintenance

1. Deflection plates swipes and DI, EtOH ok
2. Cytometer > Cleaning modes > Prepare for aseptic sort (45).
3. Change the nozzle if required (temporary solution until new integrator purchased).
4. Clean nozzle before adjusting the sort block.
 - Optional: Clean Flow Cell 3X with BD Detergent solution and DI
 - Full fluidics shutdown with eth
 - Turn off stream, purge filters.
 - Prime tanks after refill or changing filters.

Monthly maintenance

- Change waste lid white cap if necessary.

BD FACSAriaFusion™ basic SOP

BD FACSAriaFusion™ basic SOP

EXTERNAL AIR REQUIREMENTS

Air supply 6.69 Bar (96.6 psi) regulated. The source of the compressed air must deliver clean (<5 ppm total particulate) air at stable pressures. The case air consumption, based on 100 tubes per hour, is approximately 100 L per hour.

BD FACSAriaFusion™ basic SOP

Recommended Sort Settings/Parameters for BD FACSAria Fusion

	Nozzle		
	70µm	85µm	100µm
PSI	70	45	20
Amplitude	60	32	12
Frequency	87	45	30
Drop1	150	150	150
Gap (Max)	6 (14)	7 (17)	10 (21)
SSΔGap	2-3	3	4
Events per sec	18000	18000	6000